

LavaPep Interfering Compounds

At concentrations in sample greater than listed, the reagents could affect the accuracy of the LavaPep Kit. Note that a 1/400 or 1/1000 dilution in water of most sample preparations would be sufficient to dilute chemicals below the upper limit. It is recommended that samples and standards are prepared in the same buffer.

| Compound | Manufacturer | Recommended limit in LavaPep kit |
|---------------------------------|--------------|----------------------------------|
| SDS | BDH | 0.05% |
| CHAPS | BioRad | 0.01% |
| NP40 | Sigma | 0.005% |
| thiourea | Sigma | 500 mM |
| urea | BDH | 500 mM |
| triton X 100 | UT Packard | 0.005% v/v |
| tween 20 | Sigma | 0.01% v/v |
| dithiothreitol | BioRad | 1.5 mM |
| tributyl phosphine | Sigma | 5 μ M |
| Methyl methanethiol sulfonate | Sigma | 1 mM |
| Triethylammonium bicarbonate | Sigma | 2.5 mM |
| Tris (2-carboxy-thyl) phosphine | Sigma | 500 μ M |
| iodoacetamide | Sigma | 50 mM |
| calcium chloride | BDH | 500 μ M |
| tris-HCl | Sigma | 500 μ M |
| NH_4CO_3 | Sigma | 500 μ M |
| HCl | Univar | 500 μ M |
| TFA | Aldrich | 0.005% |
| Formic acid | Univar | 0.01% |
| ACN | Unichrom | 0.5% |

Table 1. Acceptable maximum levels of various chemicals and reagents for LavaPep Peptide Quantification kit.